



Auto Ref/
Keratometers

Digital
Refractors

Exam Lanes
& Furniture

Slit Lamps

Lensmeters

Microscope

AcuityCharts

Retinal Cameras

INTRODUCING **HLM-1** THE NEW LENSMETER



**FASTER
AND MORE
ACCURATE
RESULTS**



HLM-1 LENSMETER



Hartmann Sensor Wavefront Analysis Technology with more measurement points than previous generations.



Easily measures darkly-tinted or mirrored sunglasses by automatically amplifying the amount of light and calculating its refractive power without any manual configurations.



Expanded prism measurement range up to 20Δ measuring from all directions of BU, BD, BI, BO.

Wide range for measuring small or large (blank) lenses from diameters of $\varnothing 15\text{mm}$ to $\varnothing 120\text{mm}$.

Meets or exceeds the Medical Equipment Certificate, Class B standard.

Slim and compact design makes it perfect for small spaces.

Auto lens recognition including single vision and progressive lenses.

Improved progressive (multifocal) lens measurement makes measuring fast and easy by simply moving the target and following the guides on the screen.

Newly designed nose cone and lens support allows the user to measure smaller frame styles with ease.



DIAGNOSTICS



SURFACING



COATING



FINISHING



CONSUMABLES



SERVICE

