

Auto Ref/  
Keratometers

Digital  
Refractors

Exam Lanes  
& Furniture

Slit Lamps

Lensmeters

Acuity  
Charts

Ocular Surface  
Analyzers

INTRODUCING THE NEW

## IDRA OCULAR SURFACE ANALYZER



**FULLY AUTOMATED,  
COMPREHENSIVE DRY EYE  
DIAGNOSTIC SYSTEM**

**4 NON-INVASIVE TESTS  
IN 5 MINUTES**




### IDRA

OCULAR SURFACE ANALYZER




#### AUTO-INTERFEROMETRY

Auto detection and analysis providing maximum and average lipid layer thickness. Auto blinking rate.



#### TEAR MENISCUS HEIGHT MEASUREMENT

Non-invasive and quick.




#### AUTO-NIBUT (NON-INVASIVE BREAK UP TIME)

Only requires a short video. Blink rate analysis.



#### MEIBOGRAPHY AND 3D MEIBOMIAN GLAND IMAGING

Auto length and width with loss area percentage. 3D images.




VALUES ARE DISPLAYED IN A USER FRIENDLY GRADING SCALE IDEAL FOR PATIENT ENGAGEMENT



DIAGNOSTICS



SURFACING



COATING



FINISHING



CONSUMABLES



SERVICE

